

OR8B2/OR8B3 Antibody



PACO03849

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR8B2/OR8B3

Uniprot

Q96RD0/Q8NGG8

Synonyms:

OR8B2; Olfactory receptor 8B2; Olfactory receptor OR11-309; OR8B3; Olfactory receptor 8B3; Olfactory receptor OR11-311

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 8B2/3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A