

olfr89 Antibody



PACO03856

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

olfr89

Uniprot

O95499

Synonyms:

Olfactory receptor 89

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 89.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A