

# ATP5L2 Antibody



PACO03857

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

ATP5L2

**Uniprot**

Q7Z4Y8

**Synonyms:**

ATP5L2; ATP5K2; ATP synthase subunit g 2; mitochondrial; ATPase subunit g 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human ATP5L2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A