

B4GALT1 Antibody



PACO03859

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

B4GALT1

Uniprot

P15291

Synonyms:

B4GALT1; GGTB2; Beta-1; 4-galactosyltransferase 1; Beta-1,4-GalTase 1; Beta4Gal-T1; b4Gal-T1; UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 1; UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 1

Immunogen:

Synthesized peptide derived from the C-terminal region of human Beta-1,4-Gal-T1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A