

INPP5J Antibody



PACO03882

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300, IHC:1:200-1:1000

Protein Background:

Gene ID:

INPP5J

Uniprot

Q15735

Synonyms:

INPP5J; PIB5PA; PIPP; Phosphatidylinositol 4; 5-bisphosphate 5-phosphatase A; Inositol polyphosphate 5-phosphatase J

Immunogen:

Synthesized peptide derived from the C-terminal region of human INPP5J.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A