

CNTN6 Antibody



PACO03891

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CNTN6

Uniprot

Q9UQ52

Synonyms:

CNTN6; Contactin-6; Neural recognition molecule NB-3; hNB-3

Immunogen:

Synthesized peptide derived from the Internal region of human Contactin 6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A