

GNAI3 Antibody



PACO03906

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

GNAI3

Uniprot

P08754

Synonyms:

GNAI3; Guanine nucleotide-binding protein G(k) subunit alpha; G(i) alpha-3

Immunogen:

Synthesized peptide derived from the Internal region of human Galpha i-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A