

# AQP11 Antibody



PACO03914

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

AQP11

**Uniprot**

Q8NBQ7

**Synonyms:**

AQP11; AQPX1; PSEC0027; Aquaporin-11; AQP-11

**Immunogen:**

Synthesized peptide derived from the Internal region of human AQP11.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A