

STK32C Antibody



PACO03915

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

STK32C

Uniprot

Q86UX6

Synonyms:

STK32C; Serine/threonine-protein kinase 32C; PKE; YANK3

Immunogen:

Synthesized peptide derived from the Internal region of human STK32C.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A