

GPR133 Antibody



PACO03919

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

ADGRD1

Uniprot

Q6QNK2

Synonyms:

GPR133; PGR25; Probable G-protein coupled receptor 133; G-protein coupled receptor PGR25

Immunogen:

Synthesized peptide derived from the Internal region of human GPR133.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A