

OR4N4 Antibody



PACO03924

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR4N4

Uniprot

Q8N0Y3

Synonyms:

OR4N4; Olfactory receptor 4N4; Olfactory receptor OR15-1; Olfactory receptor OR15-5

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 4N4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A