

NBPF5P Antibody



PACO03936

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

NBPF5P

Uniprot

Q86XG9

Synonyms:

NBPF5; Neuroblastoma breakpoint family member 5

Immunogen:

Synthesized peptide derived from the C-terminal region of human NBPF5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A