

CYB561D1 Antibody



PACO03937

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

CYB561D1

Uniprot

Q8N8Q1

Synonyms:

CYB561D1; Cytochrome b561 domain-containing protein 1

Immunogen:

Synthesized peptide derived from the C-terminal region of human CYB561D1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A