

# RNF149 Antibody



PACO03940

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

RNF149

**Uniprot**

Q8NC42

**Synonyms:**

RNF149; DNATP2; E3 ubiquitin-protein ligase RNF149; DNA polymerase-transactivated protein 2; RING finger protein 149

**Immunogen:**

Synthesized peptide derived from the Internal region of human RNF149.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A