

GPR27 Antibody



PACO03941

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

GPR27

Uniprot

Q9NS67

Synonyms:

GPR27; SREB1; Probable G-protein coupled receptor 27; Super conserved receptor expressed in brain 1

Immunogen:

Synthesized peptide derived from the Internal region of human GPR27.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A