

GPRIN3 Antibody



PACO03947

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

GPRIN3

Uniprot

Q6ZVF9

Synonyms:

GPRIN3; KIAA2027; G protein-regulated inducer of neurite outgrowth 3; GRIN3

Immunogen:

Synthesized peptide derived from the Internal region of human GRIN3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A