

RPL7L1 Antibody



PACO03950

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

RPL7L1

Uniprot

Q6DKI1

Synonyms:

RPL7L1; 60S ribosomal protein L7-like 1

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L7L1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A