

OR13D1 Antibody



PACO03953

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR13D1

Uniprot

Q8NGV5

Synonyms:

OR13D1; Olfactory receptor 13D1; Olfactory receptor OR9-15

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 13D1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A