

FFAR2 Antibody



PACO03956

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

FFAR2

Uniprot

O15552

Synonyms:

FFAR2; GPR43; Free fatty acid receptor 2; G-protein coupled receptor 43

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR43.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A