

# TMEM30C Antibody



PACO03961

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

TMEM30CP

**Uniprot**

A0ZSE6

**Synonyms:**

TMEM30C; CDC50C; Cell cycle control protein 50C; Transmembrane protein 30C

**Immunogen:**

Synthesized peptide derived from the Internal region of human Cdc50C.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A