

# GTF3A Antibody



PACO04006

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

GTF3A

**Uniprot**

Q92664

**Synonyms:**

GTF3A; Transcription factor IIIA; TFIIIA

**Immunogen:**

Synthesized peptide derived from the Internal region of human TFIIIA.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A