

GUCY1B2 Antibody



PACO04007

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

GUCY1B2

Uniprot

O75343

Synonyms:

GUCY1B2; Guanylate cyclase soluble subunit beta-2; GCS-beta-2

Immunogen:

Synthesized peptide derived from the N-terminal region of human GCS-Beta-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A