

UBN1 Antibody



PACO04009

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

UBN1

Uniprot

Q9NPG3

Synonyms:

UBN1; Ubinuclein-1; HIRA-binding protein; Protein VT4; Ubiquitously expressed nuclear protein

Immunogen:

Synthesized peptide derived from the Internal region of human UBN-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A