

GPR171 Antibody



PACO04013

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

GPR171

Uniprot

O14626

Synonyms:

GPR171; H963; Probable G-protein coupled receptor 171; G-protein coupled receptor H963

Immunogen:

Synthesized peptide derived from the Internal region of human GPR171.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A