

GPR132 Antibody



PACO04015

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

GPR132

Uniprot

Q9UNW8

Synonyms:

GPR132; G2A; Probable G-protein coupled receptor 132; G2 accumulation protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR132.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A