

# Acetyl-H2AFZ (K5) Antibody



PACO04024

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

H2AFZ

**Uniprot**

P0C0S5

**Synonyms:**

H2AFZ; H2AZ; Histone H2A. Z; H2A/z

**Immunogen:**

Synthesized peptide derived from human Histone H2A around the acetylation site of K5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A