

KCNIP1 Antibody



PACO04044

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

KCNIP1

Uniprot

Q9NZI2

Synonyms:

KCNIP1; KCHIP1; VABP; Kv channel-interacting protein 1; KCHIP1; A-type potassium channel modulatory protein 1; Potassium channel-interacting protein 1; Vesicle APC-binding protein

Immunogen:

Synthesized peptide derived from the N-terminal region of human KCHIP1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A