

VN1R3 Antibody



PACO04055

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

VN1R3

Uniprot

Q9BXE9

Synonyms:

VN1R3; V1RL3; FKSG46; Vomeronasal type-1 receptor 3; V1r-like receptor 3

Immunogen:

Synthesized peptide derived from the Internal region of human V1RL3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A