

# HOXA1 Antibody



PACO04060

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

HOXA1

**Uniprot**

P49639

**Synonyms:**

HOXA1; HOX1F; Homeobox protein Hox-A1; Homeobox protein Hox-1F

**Immunogen:**

Synthesized peptide derived from the Internal region of human Hox-A1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A