

RAB37 Antibody



PACO04073

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

RAB37

Uniprot

Q96AX2

Synonyms:

RAB37; Ras-related protein Rab-37

Immunogen:

Synthesized peptide derived from the C-terminal region of human Rab 37.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A