

AKR1CL1 Antibody



PACO04103

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

AKR1C8P

Uniprot

Q5T2L2

Synonyms:

AKR1CL1; Aldo-keto reductase family 1 member C-like protein 1

Immunogen:

Synthesized peptide derived from the Internal region of human AKR1CL1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A