

OR2AT4 Antibody



PACO04104

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR2AT4

Uniprot

A6NND4

Synonyms:

OR2AT4; Olfactory receptor 2AT4; Olfactory receptor OR11-265

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2AT4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A