

VN2R1P Antibody



PACO04116

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

VN2R1P

Uniprot

Q8NGV9

Synonyms:

Putative calcium-sensing receptor-like 1

Immunogen:

Synthesized peptide derived from the Internal region of human VN2R1P.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A