

# GPR144 Antibody



PACO04121

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

ADGRD2

**Uniprot**

Q7Z7M1

**Synonyms:**

GPR144; PGR24; Probable G-protein coupled receptor 144; G-protein coupled receptor PGR24

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human GPR144.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A