

# DGAT2L6 Antibody



PACO04122

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

DGAT2L6

**Uniprot**

Q6ZPD8

**Synonyms:**

DGAT2L6; DC3; Diacylglycerol O-acyltransferase 2-like protein 6; Diacylglycerol O-acyltransferase candidate 3; hDC3

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human DGAT2L6.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A