

RAB41 Antibody



PACO04123

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

RAB41

Uniprot

Q5JT25

Synonyms:

RAB41; Ras-related protein Rab-41

Immunogen:

Synthesized peptide derived from the Internal region of human Rab 41.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A