

# GPR142 Antibody



PACO04133

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

GPR142

**Uniprot**

Q7Z601

**Synonyms:**

GPR142; PGR2; Probable G-protein coupled receptor 142; G-protein coupled receptor PGR2

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human GPR142.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A