

IL1A Antibody



PACO04140

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

IL1A

Uniprot

P01583

Synonyms:

IL1A; IL1F1; Interleukin-1 alpha; IL-1 alpha; Hematopoietin-1

Immunogen:

Synthesized peptide derived from the N-terminal region of human IL-1alpha.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A