

GPR153 Antibody



PACO04180

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

GPR153

Uniprot

Q6NV75

Synonyms:

GPR153; PGR1; Probable G-protein coupled receptor 153; G-protein coupled receptor PGR1

Immunogen:

Synthesized peptide derived from the Internal region of human GPR153.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A