

# OR6B2 Antibody



PACO04185

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR6B2

**Uniprot**

Q6IFH4

**Synonyms:**

OR6B2; OR6B2P; Olfactory receptor 6B2; Olfactory receptor OR2-1

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 6B2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A