

OR51Q1 Antibody



PACO04188

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR51Q1

Uniprot

Q8NH59

Synonyms:

OR51Q1; Olfactory receptor 51Q1

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 51Q1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A