

# OR5T3 Antibody



PACO04194

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR5T3

**Uniprot**

Q8NKG3

**Synonyms:**

OR5T3; Olfactory receptor 5T3; Olfactory receptor OR11-178

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 5T3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A