

# OR6C70 Antibody



PACO04201

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

OR6C70

**Uniprot**

A6NIJ9

**Synonyms:**

OR6C70; Olfactory receptor 6C70

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 6C70.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A