

OR13C2/OR13C9 Antibody



PACO04219

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR13C2/OR13C9

Uniprot

Q8NGS9/Q8NGT0

Synonyms:

OR13C2; Olfactory receptor 13C2; Olfactory receptor OR9-12; OR13C9; Olfactory receptor 13C9; Olfactory receptor OR9-13

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor O13C2/9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A