

# LGALS2 Antibody



PACO04231

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

LGALS2

**Uniprot**

P05162

**Synonyms:**

LGALS2; Galectin-2; Gal-2; Beta-galactoside-binding lectin L-14-II; HL14; Lactose-binding lectin 2; S-Lac lectin 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human Galectin-2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A