

OR2T29/OR2T5 Antibody



PACO04247

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2T29/OR2T5

Uniprot

Q8NH02/Q6IEZ7

Synonyms:

OR2T29; Olfactory receptor 2T29; OR2T5; Olfactory receptor 2T5; Olfactory receptor OR1-62

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 2T5/2T29.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A