

OR4A4P/OR4A47 Antibody



PACO04248

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR4A4P/OR4A47

Uniprot

Q8NGN8/Q6IF82

Synonyms:

OR4A4P; Putative olfactory receptor 4A4; Olfactory receptor OR11-107; OR4A47; Olfactory receptor 4A47; Olfactory receptor OR11-113

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4A4/47.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A