

MCM4 Antibody



PACO04307

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MCM4

Uniprot

P33991

Synonyms:

MCM4; CDC21; DNA replication licensing factor MCM4; CDC21 homolog; P1-CDC21

Immunogen:

Synthesized peptide derived from human MCM4 around the non-phosphorylation site of S54.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A