

**Phospho-MAP3K9/MAP3K10 (T312/266)
Antibody
PACO04335**



Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

MAP3K9/MAP3K10

Uniprot

P80192/Q02779

Synonyms:

MAP3K9; MLK1; PRKE1; Mitogen-activated protein kinase kinase kinase 9; Mixed lineage kinase 1; MAP3K10; MLK2; MST; Mitogen-activated protein kinase kinase kinase 10; Mixed lineage kinase 2; Protein kinase MST

Immunogen:

Synthesized peptide derived from human MLK1/2 around the phosphorylation site of T312/266.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A