

OR2J3 Antibody



PACO04355

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR2J3

Uniprot

O76001

Synonyms:

OR2J3; Olfactory receptor 2J3; Hs6M1-3; Olfactory receptor OR6-16; OR6-6; Olfactory receptor 6-6

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2J3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A